


# **MED-JET<sup>®</sup> - H-II**

## **PATENTED LOW PRESSURE DUAL FUNCTION INJECTOR FOR PHYSIATRICALS**

- A unique state-of-the-art needle-free system designed to function with or without a needle
- Safe, fast and extremely accurate
- Suitable for single or multiple intramuscular, intradermal and subcutaneous injections
- Adjustable and accurate dose of 0.01 to 1cc
- Propelled by CO2 or compressed air
- Adjustable air pressure to ensure patient comfort
- Can inject up to 600 injections/hour
- Saves time for the physician
- Designed for the patient's comfort: patients can be anesthetized first with the MED-JET<sup>®</sup> without the needle and then injected with the treatment drug with the needle. Thus reducing the pain associated with the needle
- Ideal for administering Cortisone<sup>®</sup>, and other pain relief injectables such as BOTOX<sup>®</sup>, Xylocaine<sup>®</sup>, etc. to patients suffering from arthritis, fibromyalgia, rheumatoid arthritis and other chronic pains
- The MED-JET<sup>®</sup> – H-II can be used with an attached needle or as thin as 30 gauge
- The advantage of using a needle with the low pressure MED-JET<sup>®</sup> injector is the fact that the drug is dispersed evenly over the affected area and thus ensuring optimization of its effect



Safe to use on most of the body parts



**THE ART AND SCIENCE OF INJECTING REINVENTED**



**MEDICAL INTERNATIONAL TECHNOLOGIES (MIT Canada) INC.**  
1872 Beaulac, Saint-Laurent, QC, Canada H4R 2E7  
Tel.: (514) 339-9355 Fax: (514) 339-2885 e-mail: [marketing@mitcanada.ca](mailto:marketing@mitcanada.ca)  
[www.mitcanada.ca](http://www.mitcanada.ca)